DentalCAM & DentalCNC 7+8 installieren (English instructions below)

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Empfehlung:

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System requirements

Recommendation:

- CPU:
- Single machine control: 64-bit Intel Core i CPU of the 10th generation with 4 cores and hyperthreading (e.g., Intel Core i5-10500)

Alternative: 64-bit AMD Ryzen CPU with at least 4 cores (e.g., AMD Ryzen 5 3600)

- Multi-machine control: 64-bit Intel Core i CPU of the 10th generation 6 cores and hyperthreading (e.g., Intel Core i7-10700)

Alternative: 64-bit AMD Ryzen CPU with at least 6 cores and 12 threads (e.g., AMD Ryzen 5 3600)

- Operating system: 64-bit version of Windows® 10, 11
- Storage device: SSD with 512 GiB
- Random access memory:
- Single machine control: 8 GiB
- Multi-machine control: 16 GiB
- Display resolution: 1920 x 1080 (Full HD)
- A permanent Internet connection during installation and use

Supported machines

dentalcam & dentalcnc support the following machines:

- E3 (only dentalcnc – use as CAM program trimcam)
- E4, E5, K4 edition, K5, K5+, N4, N4+, R5, S1, S2, S5
- Z4 with external CAM computer
- Wet grinding module NO-1

Preparing the installation

1. Install all Windows updates onto the CAM computer.
2. Update all hardware drivers.
3. If the Windows system language is neither English nor German, change the unicode and non-unicode settings to English.

Changing the unicode setting to English

Changing the non-unicode setting to English

4. If the screen resolution of the CAM computer is 1920 x 1080 (Full HD) or similar, set the text and app scaling to 100 %.

Otherwise, the application windows might not fit onto the screen.

5. Make sure that the necessary .NET Framework versions are activated:

- a. Navigate to **Control Panel > Programs & Features**.

- b. Select **Turn Windows Features on or off**.

- c. Activate **.NET Framework 3.5** and **.NET Framework 4.7**.

- d. Select **OK**.

6. Deactivate *all* energy saver settings on your computer. Keep all energy saver settings disabled while working with dentalcam and dentalcnc because otherwise the job calculation and/or machining can be aborted.

7. Disable the USB settings for selective power saving. USB selective suspend setting

8. When using a portable computer, ensure a constant electrical supply.

9. If you want to use 1 designated computer for dentalcam which will send job data to computers which are connected to the machines over your network, you may want to read the following information first:

Installing dentalcam and dentalcnc

NOTICE

Incorrect job execution and damaging of the CNC machine if the program version is outdated and/or incompatible

If you use dentalcam and/or dentalcnc versions that are outdated or incompatible with your CNC machine, jobs may be calculated or machined incorrectly or the machine may get damaged.

- Only install the latest dentalcam and dentalcnc version that is available for your CNC machine.
- Always update dentalcam and dentalcnc to the latest version that is available for your CNC machine.

1. If you received the installation data as a ZIP file, unpack it onto the hard drive of your CAM computer.
2. Disable all anti-malware software for the duration of the installation.

3. Right-click on **Setup.exe** in the unzipped folder. From the context-menu, select **Run as Administrator**.
4. If required, enter the login credentials of the administrator account.
 - The setup window opens.*The setup window: Configure a new installation*
5. Activate the **New** option.
 - At the bottom of the setup window, 2 lists are shown with available languages and supported machines.
6. If you only need dentalcnc on the computer, deactivate the **dentalcam** option.
7. Activate the **Create Desktop Icon** option. Otherwise, you will only be able to start the applications through the executable files (.exe files) in the installation folder.
8. Activate the **Install USB/Ethernet driver** option. Without this driver, the connection to your machine will not work.
9. (Optional) If you want to view STL files on your CAM computer, activate the **Install STL viewer** option.
10. (Optional) If you want to change the standard installation folder, enter the desired folder name in the **Destination folder** input field.
 - Select a folder outside of the **Program files** or **Program files(x86)** folder or otherwise the applications will not work properly.
 - Start the path with the drive letters and a colon.

The path of the installation folder may only contain the following characters:

- A–Z
- a–z
- 0–9
- _-+%@!#{ } ~&=^`
- Space

Please also note this restriction if you change the installation path later.

11. In the **Language** list, select your desired user interface language. You can change this setting in the application settings at any time.
12. In the **Install for the following machine type** list, select the machine type that you want to use the applications with.
13. To start the installation, select the following icon:
 - The setup program installs the core applications and the setup window closes. Depending on the selected components, multiple installation wizards open.
14. Install all components with the help of the installation wizards:
 - Select **Next** and **Finish** until the relevant wizard closes.
 - The wizard for the Runtime Environment for C++ may display a message that a newer version is already installed and abort the installation. Ignore this behavior, dentalcam and dentalcnc should run just fine.
15. If you install the applications for a machine with an Ethernet port, additional steps are necessary. You will find further information on this [here](#).
16. If you are operating a machine that is supported by dentalcnc 8.2, install this version. To do this, start **Setup.exe** in the installation package.
17. Activate all your anti malware applications.